



BRASSELER
USA®

XP-3D™

MINIMALLY INVASIVE ANATOMICAL CLEANING




BY YOUR SIDE

THE CHALLENGE...

NiTi instrumentation has made a significant impact on the field of endodontics. However, despite the benefit of improved efficiency, all traditional NiTi file systems have the same universal disadvantages:

► Non-anatomical shaping

The root canal system is highly complex. The majority of canals have an irregular anatomic shape. Despite significant advancements, conventional files can only make round shapes and cannot reach significant parts of the canal during treatment.




► Excessive torque/fatigue

Traditional NiTi files are prone to file deformation and potential separation, and they apply unwanted stress to the tooth.



► Require the use of multiple NiTi files

In order to reduce cyclic and torsional fatigue, multi-file sequencing techniques are utilized. This adds cost, complexity and chair-time.



Tooth image credit: *Incidence of Dentinal Defects after Root Canal Preparation: Reciprocating Versus Rotary Instrumentation. JOE. 2013; 37(4):501-504*

THE SOLUTION...


XP-3D SHAPER+

MINIMALLY INVASIVE ANATOMICAL CLEANING

The XP-3D Shaper+ is the next generation of our adaptive XP-3D line, created for shorter procedure times and agility. It addresses the shortcomings of traditional NiTi instrumentation and ushers in a new biologic standard of care in endodontic instrumentation. The XP-3D Shaper+ is a totally unique system that offers predictability, superior resistance and the ability to use multiple times in initial treatment, giving it the added advantage of shorter procedure times.^{1,2}

The patented XP-3D Shaper+ safely, efficiently and effectively cleans the root canal system three dimensionally while respecting the canal anatomy. As it rotates, the instrument's orbit expands and contracts to abrade the broad and narrow aspects of the canal equally. This intuitive micro mechanical debridement allows the practitioner to utilize a single instrument to safely and efficiently clean and enlarge the canal while respecting the original canal morphology.


Anatomically Prepared Canal




► SIZE & EXPANSION CAPACITY

The XP-3D Shaper+ is an ISO #30 with a 1° taper. Its adaptive core design allows it to adapt to canal morphology ranging from size #30 to size #90 at a .02 to .08 taper.* The flexible core gently pulsates within the root canal as it spins, continually adapting to the canal's natural anatomy.**

Tip Size
#30-#90



Taper
.02-.08




1 Azim AA, Piasecki L, da Silva Neto UX, Cruz ATG, Azim KA. XP Shaper, A Novel Adaptive Core Rotary Instrument: Micro-computed Tomographic analysis of Its Shaping Abilities 2017;43(9):15321538.doi:10.1016/j.joen.2017.04.022

2 AlOmari T, Mustafa R, Al-Fodeh R, El-Farraj H, Khaled W, Jamleh A. Debris Extrusion Using Reciproc Blue and XP Endo Shaper Systems in Root Canal Retreatment. Int J Dent. 2021;2021:6697587. Published 2021 Mar 24.doi:10.1155/2021/6697587


*For roots above size #30, the natural apical anatomy will be maintained. Taper increases ~.02 degrees with every ~10 gentle strokes. The typical taper created with XP-3D Shaper+ is .04.

**For maximum 3-dimensional cleaning potential consider also using the XP-3D Finisher.



▶ MAXWIRE®

Featuring Brasseler's exclusive MaxWire® Technology, the XP-3D Shaper+ adapts to the canal's natural anatomy by expanding once exposed to body temperature.



MaxWire


- Super elasticity, extreme flexibility and resistance to cyclic fatigue
- Transformation to a robust, predefined serpentine shape at body temperature 95°F (35°C)
- Ability to gently expand to the natural anatomy of the root canal

68°F (20°C)

M-PHASE

Martensitic Phase

At or below room temperature the instrument is very malleable and has a relaxed serpentine shape.




Soft Shape

95°F (35°C)

A-PHASE

Austenitic Phase

When introduced to warmer temperatures (>95°F), the instrument transitions to a more robust serpentine shape.



Memorized Shape

▶ PRODUCT OFFERING (3 PACK)

File/Tip ID	Tip Size	21 mm (Stopper)	25 mm (Stopper)	31 mm (Stopper)
▲ Blue	30	5028597U0	5028598U0	5028599U0

Patent Pending


▶ BOOSTER TIP™

The new XP-3D Shaper+ is equipped with an improved ShaperTip design which helps guide the serpentine XP-3D Shaper+ around curvatures and keeps it centered in the canal. The tip of the XP-3D transitions from an ISO# 15 to an ISO# 30 within 1mm from the tip, thereby functioning as both a scouting and finishing file.

- Six facets and six cutting edges for optimal guidance and scouting
- Shaping begins at ISO #15 and transitions to ISO #30
- Tip transitions from six to three cutting edges to improve the clearing of debris


Original XP-3D Booster Tip

3 Facets



XP-3D Shaper+ Improved Tip

6 Facets




▶ ADAPTIVE CORE™ TECHNOLOGY

The XP-3D Shaper+ features unique Adaptive Core Technology which allows the smaller central core of the file to move freely and adapt to the canal's natural morphology. This facilitates debris removal, making it more efficient without occluding the dentinal tubules. Furthermore, the turbulence generated by the XP-3D Shaper+ enhances the penetration of irrigants into dentinal tubules and improves the overall disinfection of the canal.


Traditional Core Files

Debris Accumulation/Compaction



XP-3D Shaper+ Adaptive Core

Free-Floating with Turbulance




● Debris

● Stress applied to the canal wall

▶ GENTLE, LOW-STRESS MECHANICAL DEBRIDEMENT

The XP-3D Shaper+'s design characteristics drastically limit the amount of torque and stress applied to both the instrument and the canal. This results in reduced instrument separation and dentinal micro-cracks.

The below graph shows the result of a test performed on plastic blocks demonstrating the average torque generated by three instruments during canal treatment.




● XP-3D Shaper+ ● RaCe (Brasseler) ● Competitor

On average, the XP-3D Shaper+ generates 47% less torque than the RaCe instrument and 88% less than the leading competitor's instrument.*


The images below show the result of a photoelasticity test performed on plastic blocks to illustrate the stress applied by instruments during canal treatment. The greater the stress, the higher the risk of micro-cracks and instrument separation.

Competitor A*




*Alternating Movement:
Very High Stress*

Competitor B*



*Rotary Movement:
High Stress*


XP-3D™ Shaper+*




*Rotary Movement:
Very Low Stress*

*Internal testing available upon request.


▶ SUPERIOR RESISTANCE*



▶ SUPERIOR FLEXIBILITY*



▶ HIGH BACTERIA REDUCTION - UP TO 99.5%³



*Internal testing available upon request.


³ Siddique R, Nivedhitha MS, Ranjan M, Jacob B, Solete P. Comparison of antibacterial effectiveness of three rotary file system with different geometry in infected root canals before and after instrumentation-a double-blinded randomized controlled clinical trial. *BDJ Open*. 2020;6:8. Published 2020 Jun 8. doi:10.1038/s41405-020-0035-7. <http://creativecommons.org/licenses/by/4.0/>. No changes made.

XP-3DTM

SCOUT FILE

The XP-3D Scout File is the newest offering in Brasseler's innovative XP-3D file system. With its small size and precise design, it offers clinicians an improved ability to scout and prepare root canals during advanced cases.

Available in sizes #15/.04 in 21 and 25mm lengths, this file is designed to be used in conjunction with the new XP-3D Shaper+ and XP-3D Finisher files during root canal therapy. By incorporating the XP-3D Scout File during your more advanced cases, you can more effectively navigate through curved and narrow canals, providing increased efficiency and accuracy during the procedure.



▶ PRODUCT OFFERING (3 PACK)

File/Tip ID	Tip Size	21 mm (Stopper)	25 mm (Stopper)
△ White	15	5028594U0	5028593U0

XP-3DTM


PROKIT

For simpler ordering and convenience, try XP-3D ProKit, available in 21 and 25mm. Each kit contains everything needed for your full XP-3D root canal therapy.

- Kit Contents:
- K-File #10/.02
 - XP-3D Scout #15
 - XP-3D Shaper+ #30
 - XP-3D Finisher #30

▶ PRODUCT OFFERING (4 PACK)

File/Tip ID	Tip Size	21 mm (Stopper)	25 mm (Stopper)
Asst.	Asst.	5028596U0	5028595U0




▶ CLINICAL APPLICATIONS

- Anatomical debridement of basic and advanced root canal anatomy
- Retreatment of failed root canals
- Enhanced irrigation


▶ CLINICAL TIPS

- Packaged sterile
- 800-1000 RPMs; Torque: 1Ncm
- Use long gentle strokes to progress down to working length (WL). If WL is not reached in 8 strokes, stop, irrigate and recapitulate.
- Irrigate frequently

BASIC CASES



ADVANCED CASES



Final Preparation Size: Once WL is reached, a minimal preparation of #30/.02 is achieved (~5 strokes). The taper can be increased to .04 with an additional 10 strokes. The tip will not expand beyond #30 unless the apex is naturally larger than #30. The XP-3D Shaper+ tip can reach up to size #90 if space permits.

Prior to obturation, it is recommended to confirm the final apical diameter with an apical verifier.

- Recommended: XP-3D Apical Verifier, now available in sizes #25-50/.04 and in assortment packs.

Obturate with bioceramic bonded obturation (BC Sealer/Filler™ with BC Points™ or normal gutta percha).

XP-3DTM

APICAL VERIFIER

Available in sizes #25-50 in .04/25mm length, XP-3D Apical Verifiers allow you to easily confirm the final apical diameter and take any guesswork out of obturation.

- Sterilizable
- Ni-Ti construction
- Color-coded for easy identification
- Available in sizes #25-50 in .04 taper 25mm length as well as an assortment pack of 25-50




▶ PRODUCT OFFERING (6 PACK)

File/Tip ID	Tip Size	25 mm (Stopper)
▲ Red	25	5030079U0
▲ Blue	30	5030080U0
▲ Green	35	5030081U0
▲ Black	40	5030082U0
△ White	45	5030083U0
▲ Yellow	50	5030084U0
Assortment	25-50	5030085U0

CASE STUDIES

► CASE 1

Canal preparation (ex-vivo) of a maxillary right first premolar to size #30/.04 with the XP-3D Shaper then filled with EndoSequence® BC Sealer™ and EndoSequence® BC Points™. We can see that the original shape of the canal has been perfectly preserved.




Radiograph showing the bucco-lingual aspect of the maxillary first premolar



Cross-section 1mm from the apex



Cross-section 4mm from the apex



Cross-section 7mm from the apex

▶ CASE 2

A 62-year-old woman presenting with symptomatic pulpitis on the upper right first molar. After preparing a glide path to #15/.02, the canals were instrumented with the XP-3D Shaper to #30/.04 following the instructions for use. The canals were then obturated with EndoSequence® BC sealer™ and EndoSequence® BC Points™.



Pre-op



Post-op


© Dr. Kleber K.T. Carvalho (Brazil). All rights reserved.

▶ CASE 3


Pulpectomy performed on lower right first molar. After preparing a glide path, the five canals were instrumented with the XP-3D Shaper. The canals were further cleaned with the XP-3D Finisher and obturated with EndoSequence® BC Sealer™ and EndoSequence® BC Points™.



Pre-op



Post-Op: View of canals after instrumentation with the XP-3D Shaper and bioceramic obturation




Microscopic view (x12) of 3 mesial canals after instrumentation with the XP-3D Shaper, and cleaning with the XP-3D Finisher




Microscopic view (x12) of 3 mesial canals after obturation with EndoSequence BC Sealer and EndoSequence BC Points.

XP-3DTM FINISHER

MINIMALLY INVASIVE
ANATOMICAL CLEANING



Like the XP-3D Shaper+, the XP-3D Finisher utilizes Brasseler's exclusive MaxWire Technology to adapt to the canal's natural anatomy. The XP-3D Finisher has a bowed shape at body temperature and is incredibly flexible. The instrument debrides the root canal system 3-dimensionally and allows for enhanced irrigation. Unlike the XP-3D Shaper+, the XP-3D Finisher is intended to clean a prepared canal and will not change the shape once prepared.



20°C
M-PHASE
Martensitic Phase

Soft Shape



35°C
A-PHASE
Austenitic Phase

Memorized Shape




▶ PRODUCT OFFERING (3 PACK)

File/Tip ID	Tip Size	21 mm (Stopper)	25 mm (Stopper)
▲ Red	25	5025842U0	5024933U0
▲ Blue	30	5025844U0	5025843U0


Patent Pending

► SIZE & EXPANSION CAPACITY

The XP-3D Finisher is an ISO #25 with a 0° taper. Its capacity to expand improves its reach 100-fold compared to a standard instrument.



Standard Instrument



XP-3D Finisher

► CLINICAL APPLICATIONS

- 3-D debridement
- Enhanced irrigation
- Retreatment

► CLINICAL TIPS


- 800-1000 RPMs; Torque:1Ncm
- XP-3D Finisher should be used only after canal preparation to at least #25
- In multirooted teeth, start with the largest canal
- Work along the entire length of the canal for approximately one (1) minute
- Place the XP-3D Finisher into the orifice of the canal prior to placing the irrigant and begin rotation

► PACKAGING

Three instruments in a sterile blister pack, for single use (each instrument can be used to clean one tooth with up to four canals).




The instruments are stored inside a plastic tube so their straight shape can be maintained or restored and the working length can be defined.



IN VITRO REVIEW


▶ BENCH-TOP CROSS-SECTION REVIEW

Mesial root of an artificial lower molar sectioned at 1, 3, 5 and 7mm from apex:




▶ BENCH-TOP CLEANING & OBTURATION SIMULATION


1 Pre-op with pulp.




2 During preparation with XP-3D Finisher.



3 After preparation with XP-3D Finisher.



4 After obturation with EndoSequence® BC Sealer™ and EndoSequence® BC Points™.





CLINICAL CASES

▶ CASE 1


Micro CT of the Distal and Mesial roots of a lower molar instrumented to 35/.04 with round NiTi files and then after final cleaning with the XP-3D Finisher.

- Pre-Op and Post-Op NiTi Pictures: Show debris in the canal and in the isthmus areas.
- Post-Op XP-3D Finisher Pictures: After final cleaning with the XP-3D Finisher, no debris is seen.



▶ CASE 2

Internal resorption case demonstrating the expansive potential of the XP-3D Finisher. Cases with irregular shaped canals and roots which are structurally compromised are ideal beneficiaries of the new XP-3D Finisher. This new technology allows us to preserve critical root structure while thoroughly disinfecting the canal.




BIOCERAMIC BONDED OBTURATION

In the era of modern endodontics, shaping is no longer dictated by the limitations of obturation materials. The introduction of non-shrinking bonded obturation (BC Sealer™ and BC Points™) allows practitioners to embrace the adaptive and minimally invasive design of the XP-3D Shaper+. Unlike traditional sealers, BC Sealer does not shrink and therefore it is not necessary to condense gutta percha in an effort to minimize the sealer interface. With BC Sealer the function of gutta percha is simply to take up space, provide a path for retreatment and provide for hydraulics/delivery of the sealer.

BC Sealer/Filler™

- Highly Biocompatible
- Absolutely Zero Shrinkage
- Produces Hydroxyapatite
- Highly Antibacterial (+12pH)
- Bonds to Dentin and BC Points™



▶ XP-3D™ OBTURATION TECHNIQUE

1. Confirm the apical diameter with XP-3D apical verifier, gutta percha point, paper point or file.
2. Coat the canal walls with BC Sealer.
3. Coat the appropriate gutta percha point with BC Sealer and place to length.
4. Sear off gutta percha at the canal orifice and vertically compact with the appropriately sized plugger.

- ▶ 1. Bao P, Shen Y, Lin J, Haapasalo M. In Vitro Efficacy of XP-3D Finisher with Two Different Protocols on Biofilm Removal from Apical Root Canals. *JOE*. 2016; Dec 7. DOI: 10.1016/j.joen.2016.09.021.

Subject: The purpose of this study was to evaluate the effectiveness of the XP-3D Finisher in biofilm removal in comparison with conventional needle irrigation (CNI) and passive ultrasonic irrigation (PUI) using an infected tooth model with an artificial apical groove.

Results/Conclusions: XP-3D Finisher showed the best biofilm removal efficacy inside and outside the groove followed by PUI and CNI ($P < .05$).

Significance: Biofilm removal is one of the most challenging phases in endodontic instrumentation. In this study the XP-3D was most effective compared to the most popular methods used to this point.

- ▶ 2. Keskin C, Sariyilmaz E, Sariyilmaz O. Efficacy of XP-3D Finisher File in Removing Calcium Hydroxide from Simulated Internal Resorption Cavity. *JOE*. 2017; Jan. DOI: 10.1016/j.joen.2016.09.009.

Subject: The aim of this study was to evaluate the effect of supplementary use of XP-3D Finisher file, passive ultrasonic activation (PUI), EndoActivator (EA), and CanalBrush (CB) on the removal of calcium hydroxide (CH) paste from simulated internal resorption cavities.

Results/Conclusions: XP-3D Finisher and PUI removed significantly more CH than SI, EA, and CB ($P < .05$), showing no significant difference between them ($P > .05$). Differences among SI, EA, and CB were also non-significant ($P > .05$).

Significance: Removal of calcium hydroxide is an indirect indicator of cleaning effectiveness against microbes. XP-3D Finisher was equal to PUI and superior to other traditional methods.

There was no indication in this study that the teeth or irrigating solution were heated above body temperature, which is required for optimal effectiveness of the XP-3D Finisher.

- ▶ 3. Alves FR, Andrade-Junior CV, Marceliano-Alves MF, Perez AR, Racas IN, Versiani MA, Sousa-Neto MD, Provenzano JC, Siquera JF JR. Adjunctive Steps for Disinfection of the Mandibular Molar Root Canal System: A Correlative Bacteriologic, Micro-Computed Tomography, and Cryopulverization Approach. *JOE*. 2016; Nov. DOI: 10.1016/j.joen.2016.08.003

Subject: This study evaluated the disinfecting ability of chemomechanical preparation with rotary nickel-titanium instruments, followed by 2 distinct adjunctive procedures in the root canals of extracted mandibular molars by means of a correlative analytical approach.

Results/Conclusions: Both XP-3D Finisher and passive ultrasonic irrigation exhibited antibacterial effectiveness, but only the former caused a significant reduction in the bacterial counts after chemomechanical preparation.

Significance: Showed superior antibacterial effectiveness compared to PUI which has been considered “state of the art”.

- ▶ 4. Wigler R, Dvir R, Weisman A, Malalon S, Kfir A. Efficacy of XP-3D Finisher files in the removal of calcium hydroxide paste from artificial standardized grooves in the apical third of oval root canals. *International Endodontic Journal*. 2016, Jun. DOI: 10.1111/iej.12668.

Subject: To compare the efficacy of the XP-3D Finisher file (XP) (FKG Dentaire, La Chaux de Fonds, Switzerland) to that of passive ultrasonic irrigation (PUI) and conventional syringe and needle irrigation (SNI) in the removal of calcium hydroxide paste from an artificial standardized groove in the apical third of root canals.

Results/Conclusions: XP-3D and PUI were more effective in removing Ca(OH)₂ from artificial standardized grooves in the apical third of root canals than SNI.

Significance: Removal of calcium hydroxide is an indirect indicator of cleaning effectiveness against microbes. XP-3D Finisher was equal to PUI and superior to other traditional methods.

There was no indication in this study that the teeth or irrigating solution were heated above body temperature, which is required for optimal effectiveness of the XP-3D Finisher.

- ▶ 5. Alves F, Marcelano-Alves M, Sousa JC, Silveira S, Prvenzano J, Siqueira J. Removal of Root Canal Filling in Curved Canals Using Either Reciprocation Single- or Rotary Multi-Instrument System and a Supplementary Step with the XP-3D Finisher. *JOE*. 2016; 42(7). DOI: 10.1016/j.joen.2016.04.007.

Subject: This study compared the efficacy of a reciprocating single-instrument system and a rotary multi-instrument system followed by a supplementary approach with a finishing instrument in removing the filling material from curved canals during retreatment.

Results/Conclusions: The rotary multiple-instrument system was more effective and faster than the reciprocating single-instrument approach in removing previous root canal fillings. The adjunctive finishing instrument XP-3D Finisher significantly improved filling material removal.

Significance: The XP-3D Finisher resulted in additional cleaning after the use of round files for retreatment purposes.

- ▶ 6. Elnaghy AM1, Mandorah A2, Elsaka SE. Effectiveness of XP-3D Finisher, EndoActivator, and File agitation on debris and smear layer removal in curved root canals: a comparative study. *Odontology*. 2016, May. DOI: 10.1007/s10266-016-0251-8.

Subject: The purpose of this study was to assess the efficacy of the XP-3D Finisher (FKG Dentaire SA, La Chaux-de-Fonds, Switzerland) file on debris and smear layer removal in curved root canals in comparison to different irrigation regimens.

Results/Conclusions: Irrigation of curved root canals using XP-3D Finisher and EndoActivator methods appears to be more effective on debris and smear layer removal than the other tested groups.

Significance: XP-3D Finisher shows value in removal of debris and smear layer, which should contribute to improved success.

There was no indication in this study that the teeth or irrigating solution were heated above body temperature, which is required for optimal effectiveness of the XP-3D Finisher.

- ▶ 7. Azim A, Aksel H, Zhuang T, Mashtare T, Babu J, Huang G. Efficacy of Four Irrigation Protocols in Killing Bacteria Colonized in Dentinal Tubules Examined by a Novel Confocal Laser Scanning Microscope Analysis. *JOE*. 2016; 42(6). DOI: 10.1016/j.joen.2016.03.009.

Subject: The aim of this study was to determine the efficiency of 4 irrigation systems in eliminating bacteria in root canals, particularly in dentinal tubules.

Results/Conclusions: XP-3D appears to be more efficient than other 3 techniques in disinfecting the main canal space and up to 50 µm deep into the dentinal tubules.


Significance: Microbial reduction is the key to endodontic success. Therefore this result shows the XP-3D superiority in that regard.

- ▶ 8. Leoni GB, Versiani MA, Silva-Sousa YT, Brniera JFB, Pecora JD, Sousa-Neto MD. Ex vivo evaluation of four final irrigation protocols on the removal of hard tissue debris from the mesial root canal system of mandibular first molars. *International Endodontic Journal*. 2016, Apr. DOI:10.1111/iej.12630.

Subject: To evaluate the efficacy of four final irrigation protocols on the reduction of hard-tissue debris accumulated within the mesial root canal system of mandibular first molars using micro-CT analysis.

Results/Conclusions: The PUI technique and XP-3D Finisher instruments were associated with significantly lower levels of AHTD compared with conventional irrigation and modified SAF system protocol in the mesial root canals of mandibular molars. The PUI method was most effective in the coronal half of the canal while the XP-3D Finisher was most effective in the apical half of the canal.

Significance: Less debris particularly in the apical half of the canal is associated with higher success.

- 
9. Sanabria-Liviach D, Moldauer BI, Garcia-Godoy F, Antonio-Campos A, Casaretto M, Torres-Navarro J, Scalerio JM. Comparison of the XP-3D Finisher File System and Passive Ultraonic Irrigation (PUI) on Smear Layer Removal after Root Canal Instrumentation Effectiveness of Two Irrigation Methods on Smear Layer Removal. *Journal of Dental and Oral Health*. 2017; 4: 1-7.


Subject: The purpose of this study was to evaluate the effectiveness of two different final irrigant activation methods in removing the smear layer at 3 and 7 mm from the apex.

Results/Conclusions: The XP-3D Finisher NaOCl/EDTA group was more effective than PUI with and without EDTA for smear layer removal at both middle and apical canal thirds.


Significance: Smear layer removal is thought to be critical for optimal success in endodontics.



DIAMONDS



CARBIDES



POLISHERS



PROCEDURE SYSTEMS



BIOCERAMICS



Since 1976, our absolute focus has been to develop products and provide services that support the practice of core operatory procedures. Today, Brasseler USA® is the premier dental instrumentation company in North America offering the most extensive selection of dental and surgical instrumentation under one brand.



HANDPIECES



ENDODONTICS



SMALL EQUIPMENT



HYGIENE



HAND INSTRUMENTS



ORAL SURGERY



LABORATORY INSTRUMENTS



Visit our website at BrasselerUSA.com

To order call 800.841.4522 or fax 888.610.1937. In Canada call 800.363.3838

